






Introducing CURVEO™

A totally new gluteus implant,
adapted from **IMPLEO™** gel.

-  A highly cohesive, form stable gel for a **natural look**.
-  New orientation marks¹ allow for **more precision during placement**.
-  Optimised ratio between projection and base width,
for **better aesthetic results**.



NAGÔR™

Reasons to choose the NEW CURVEO™ gluteal implants.

A HIGHLY COHESIVE, FORM STABLE GEL FOR A NATURAL LOOK.

Developed using gel adapted from the bestselling IMPLEO™ breast implant.

It has a thick elastomer shell that's resilient but also highly elastic, ensuring exceptional product performance.

NEW ORIENTATION MARKS¹ ALLOW FOR MORE PRECISION DURING PLACEMENT.

New orientation dots¹ and orientation line markers, for more precision during placement.

The implant benefits from an opaque, smooth shell surface to minimise abrasion.

OPTIMISED RATIO BETWEEN PROJECTION AND BASE WIDTH, FOR BETTER AESTHETIC RESULTS.

CURVEO™ has two base configurations, with a total of 13 variations – six-Round and seven-Anatomical – to provide more choice.

The optimised ratio between projection and base width creates balanced and naturally aesthetic results, when used either intramuscular or subfascially.

Round Vs Anatomical

Round Implant

An orientation dot indicates the centre of the implant and Point of Maximum Projection.




Anterior surface

Posterior surface

Anatomical Implant

Horizontal orientation line marks the centre of the implant.



Anterior surface

Posterior surface

Six orientation dots:

Three dots on Anterior surface; two to indicate inferior and one to indicate superior point of the implant. Three dots to Posterior to match Anterior surface.

GLT-R Round Gluteal Implant

Product Code	Diameter (mm) ±2	Projection (mm) ±2	Weight (g) ±1g
GLT-R 200	111	33	200
GLT-R 240	116	35	240
GLT-R 260	119	36	260
GLT-R 300	124	38	300
GLT-R 330	130	38	330
GLT-R 380	136	41	380

GLT-A Anatomical Gluteal Implant

Product Code	Height (mm) ±2	Width (mm) ±2	Projection (mm) ±2	Weight (g) ±1g
GLT-A 215	136	99	28	215
GLT-A 285	139	109	33	285
GLT-A 305	143	110	34	305
GLT-A 345	145	114	36	345
GLT-A 385	150	117	42	385
GLT-A 440	152	121	46	440
GLT-A 470	152	121	48	470

For more information about our products please visit:

www.gcaesthetics.com or email info@gcaesthetics.com

Global Consolidated Aesthetics Limited, registered in Ireland with number 450181.

Registered office: Suite 601, Q House, Furze Road, Sandyford Industrial Estate, Dublin 18, Ireland. © Copyright 2016 GC Aesthetics A50-168.

¹ Compared to other GC Aesthetics gluteus implants.