

Noa

By GC Aesthetics™

PRODUCT CATALOGUE



After-care

by GC Aesthetics

As a global market leader with 40 years of experience in medical aesthetics, GC Aesthetics is passionate about moving the industry forward through innovation and extraordinary customer experiences.

We believe in supporting our community and the patient journey.

Recovering from surgery

Clinical studies show that compression garments aid the body's healing process from procedures that cause trauma to the body.¹

We see compression as an integral part of your patients' treatment. In fact, we believe recovery after surgery is just as important as the surgery itself.

Introducing Noa

Best compression, better results.

WHAT IS NOA?

A new range of medical-grade compression garments to ensure optimum results post-surgery.

Developed in partnership with Marena, a leading garment manufacturer for 25 years, Noa's range is manufactured in the US and exceeds all industry quality standards, including FDA and CE.

Engineered for even pressure on the skin, our patented 3-dimensional stretch fabric (TriFlex™) provides optimal support and grip for a secure fit.²

Our fabric is extremely soft, and stretches 100% more than competitors, making our garments more comfortable to wear. Comfort that improves patient compliance and helps optimise post-surgical results.^{3,4}

What's more, our garments are available in a wide range of sizes that adjust to the face, breast and body contours of your patients.

FACE

NOA 901


PAGES 6-9



BREAST

NOA 003 007 510

PAGES 10-21



BODY

NOA 310 147 157

PAGES 22-33



FACE

NOA-901

Unisex Minimal
Coverage Face Mask

PROCEDURES

Oral Surgery
Maxillofacial Surgery



NOA-901

Unisex Minimal Coverage Face Mask

COLOURS ● SIZES S-XL

Clinically designed for comfort and ease of use, the face mask features adjustable hook-and-loop closures at the top and base of the head, for effortless post-operative placement.



Product Features

- Adjustable hook-and-loop closure at top of head.
- Adjustable hook-and-loop closure at base of head.

NOA-901

Unisex Minimal Coverage Face Mask

Sizing Information

SIZE	NECK SIZE
S	33-38 cm
M	38-43 cm
L	43-48 cm
XL	48-53 cm

How to Measure

Our sizing utilises more explicit criteria than normal garment sizing. Accurate sizing is critical for optimal fit and comfort. Please check your measurements and use the size chart above to determine the correct size.



NECK
Measure the circumference of your neck

Packaging Info

Weight 0.11 kg
Dimensions 17 cm x 22 cm
For product codes see pages 34-35

BREAST

NOA-003

Seamless Cup Basic Bra



PROCEDURES

Augmentation Mammoplasty
Reduction Mammoplasty
Mastopexy



NOA-003

Seamless Cup Basic Bra

COLOURS   SIZES 32-42 B/C 34-42 D/E

Clinically designed for comfort and ease of use. The adjustable hook-and-eye front closure and shoulder straps allow for effortless intraoperative placement, with a special band positioning below the inframammary fold to avoid incision interference and irritation.



Product Features

- Microfiber adjustable hook-and-eye closures on shoulder straps and front.
- Smooth, advanced seam design to remove pressure and enhance comfort.
- High back and underarm area for smooth fit and comfort
- Unique FlexFit™ cup designed to provide consistent support with changes in breast volume.
- Extended length to avoid inframammary fold.

NOA-003

Seamless Cup Basic Bra

Sizing Information

UNDER BUST	BUST	
68.5-71 cm (UK 32)	86-89 cm (B-C)	
73.5-76 cm (UK 34)	91.5-94 cm (B-C)	96.5-99 cm (D-E)
79-81 cm (UK 36)	96.5-99 cm (B-C)	102-104 cm (D-E)
84-86.5 cm (UK 38)	102-104 cm (B-C)	107-109 cm (D-E)
89-91.5 cm (UK 40)	107-109 cm (B-C)	112-114.5 cm (D-E)
94-96.5 cm (UK 42)	112-114 cm (B-C)	117-119.5 cm (D-E)

How to Measure

Our sizing utilises more explicit criteria than normal bra sizing.

Follow the steps below to measure as accurately as possible:



STEP 1

Measure your under bust to find your band size.



STEP 2

Measure the fullest part of your bust to find your cup size.

Packaging Info

Weight 0.115kg

Dimensions 17 cm x 22 cm

For product codes see pages 34-35

BREAST

NOA-007

Seamless Cup Bra
with Implant Stabiliser

PROCEDURES

Augmentation
Mammoplasty
Implant Stabilisation



NOA-007

Seamless Cup Bra with Implant Stabiliser

COLOURS   SIZES 32-42 B/C 34-42 D/E

Clinically designed for comfort and ease of use. The adjustable hook-and-eye front closure and shoulder straps allow for effortless intra-operative placement, with a special band positioning below the inframammary fold to avoid incision interference and irritation. Added stabiliser band for enhanced implant stability.



Product Features

- Microfiber 3-row hook-and-eye closures on shoulder straps.
- Built-in implant stabiliser band with multiple positioning guides.
- Soft material to avoid irritation under arm.
- Unique FlexFit™ cup designed to provide consistent support with changes in breast volume.
- Microfiber 3-row hook-and-eye front closure.
- Smooth elastic under bust band.

NOA-007

Seamless Cup Bra with Implant Stabiliser

Sizing Information

UNDER BUST	BUST	
68.5-71 cm (UK 32)	86-89 cm (B-C)	
73.5-76 cm (UK 34)	91.5-94 cm (B-C)	96.5-99 cm (D-E)
79-81 cm (UK 36)	96.5-99 cm (B-C)	102-104 cm (D-E)
84-86.5 cm (UK 38)	102-104 cm (B-C)	107-109 cm (D-E)
89-91.5 cm (UK 40)	107-109 cm (B-C)	112-114.5 cm (D-E)
94-96.5 cm (UK 42)	112-114 cm (B-C)	117-119.5 cm (D-E)

How to Measure

Our sizing utilises more explicit criteria than normal bra sizing.

Follow the steps below to measure as accurately as possible:



STEP 1


Measure your under bust to find your band size.



STEP 2

Measure the fullest part of your bust to find your cup size.

How to Use



First: Slide each arm into the arm holes or sleeves. Close the hook and eyes and/ or zipper in front, use the middle row or the most comfortable setting to start. After adjusting the front close, adjust the shoulder straps to be even.

Second: Position left side band so that the top edge of the band aligns with the top edge of the positioning dot.

Third: Position right side band so that the top edge of the band aligns with the top edge of the positioning dot, secure hook tape to left side band.

Packaging Info

Weight 0.12 kg Dimensions 17 cm x 22 cm
For product codes see pages 34-35

Manufactured by

 Marena

BREAST

NOA-510

Implant Stabiliser Band

PROCEDURES

Augmentation
Reconstruction
Mammoplasty



NOA-510

Implant Stabiliser Band

COLOURS ONE SIZE Up to 122 cm bust

This adjustable elastic implant stabiliser band is made of supportive, soft elastic with a hook-and-loop closure for use after breast augmentation procedures.



Product Features

- Adjustable hook-and-loop closure.
- Soft, elastic band.
- Fits up to a 122cm bust.

NOA-510

Implant Stabiliser Band



How to Measure

Our sizing utilises more explicit criteria than normal garment sizing. This implant stabiliser band is adjustable and fits up to a 122 centimetre bust size.

Packaging Info

Weight 0.27 kg

Dimensions 17 cm x 22 cm

For product codes see pages 34-35

BODY

NOA-310

Women's Vest
with Open Bust

PROCEDURES

Brachioplasty
Arm Contouring



NOA-310

Women's Vest with Open Bust

COLOURS   SIZES 2XS-4XL

Clinically designed for comfort and ease of use. This open-bust vest features adjustable hook-and-eye closures that allow for effortless intra-operative placement.



Product Features

- 3-row hook-and-eye upper front closure.
- Open bust.
- 3-row hook-and-eye lower front closure.
- Soft, elastic band.
- ¾-length sleeves.



NOA-310

Women's Vest with Open Bust

Sizing Information

SIZE	ARM	UNDER BUST
2XS	22-24 cm	58-64 cm
XS	24-26 cm	64-69 cm
S	26-27 cm	69-74 cm
M	27-29 cm	74-79 cm
L	29-30 cm	79-84 cm
XL	30-33 cm	84-91 cm
2XL	33-36 cm	91-102 cm
3XL	36-40 cm	102-112 cm
4XL	40-44 cm	112-122 cm

How to Measure


Our sizing utilises more explicit criteria than normal bra sizing.

Follow the steps below to measure as accurately as possible:



STEP 1

Measure the widest part of your arm.



STEP 2

Measure your under bust.

Packaging Info

Weight 0.17 kg

Dimensions 17 cm x 22 cm

For product codes see pages 34-35

BODY

NOA-147

Women's Girdle
with High-Back

CALF-LENGTH

PROCEDURES

Liposuction
Abdominoplasty



NOA-147

Women's Girdle with High-Back CALF-LENGTH

COLOURS  SIZES S-3XL

Designed for comfortable support and ease of use after liposuction and body procedures, this bodysuit is adjustable through fluctuations in swelling, and features slide-adjustable shoulder straps and graduated compression in the lower leg.



Product Features

- Adjustable shoulder straps.
- High-back.
- Elevated strap platform for extra comfort under arm.
- Adjustable 3-row hook-and-eye closure.
- Molded buttocks.
- Open crotch.
- Clean-cut finish on legs to prevent binding.
- Graduated compression in the legs.

NOA-147

Women's Girdle with High-Back CALF-LENGTH

Sizing Information

SIZE	WAIST	HIPS	THIGH
S	66-71 cm	91-97 cm	51-54 cm
M	71-76 cm	97-102 cm	54-57 cm
L	76-81 cm	102-107 cm	57-60 cm
XL	81-89 cm	107-114 cm	60-65 cm
2XL	89-99 cm	114-124 cm	65-70 cm
3XL	99-109 cm	124-135 cm	70-77 cm

How to Measure

Our sizing utilises more explicit criteria than normal garment sizing.

Follow the steps below to measure as accurately as possible. Note, the point of procedure determines the first measurement site.



WAIST

Measure the waist; identify by bending to the side to find the natural crease.



HIPS

Measure the fullest part of the hips and buttocks.



THIGH

Measure the widest point of the thigh, near the crotch.

Packaging Info

Weight 0.32 kg

Dimensions 17 cm x 22 cm

For product codes see pages 34-35

Manufactured by

 Marena

BODY

NOA-157

Women's Girdle
with High-Back

SHORT-LENGTH

PROCEDURES

Liposuction
Abdominoplasty



NOA-157

Women's Girdle with High-Back SHORT-LENGTH

COLOURS  SIZES 2XS-4XL

Designed for comfortable support and ease of use after liposuction and body procedures, this short-length bodysuit is adjustable through fluctuations in swelling, and features slide-adjustable shoulder straps.



Product Features

- Adjustable shoulder straps.
- High-back.
- Elevated strap platform for extra comfort under arm.
- Adjustable 3-row hook-and-eye closure.
- Molded buttocks.
- Open crotch.
- Clean-cut finish on legs to prevent binding.

NOA-157

Women's Girdle with High-Back SHORT-LENGTH

Sizing Information

SIZE	WAIST	HIPS	THIGH
2XS	56-61 cm	81-86 cm	44-48 cm
XS	61-66 cm	86-91 cm	48-51 cm
S	66-71 cm	91-97 cm	51-54 cm
M	71-76 cm	97-102 cm	54-57 cm
L	76-81 cm	102-107 cm	57-60 cm
XL	81-89 cm	107-114 cm	60-65 cm
2XL	89-99 cm	114-124 cm	65-70 cm
3XL	99-109 cm	124-135 cm	70-77 cm
4XL	109-119 cm	135-145 cm	77-86 cm

How to Measure

Our sizing utilises more explicit criteria than normal garment sizing.

Follow the steps below to measure as accurately as possible. Note, the point of procedure determines the first measurement site.



WAIST

Measure the waist; identify by bending to the side to find the natural crease.



HIPS

Measure the fullest part of the hips and buttocks.



THIGH

Measure the widest point of the thigh, near the crotch.

Packaging Info

Weight 0.27 kg

Dimensions 17 cm x 22 cm

For product codes see pages 34-35

Manufactured by

 Marena

Product Codes

NOA-901

Unisex Minimal Coverage Face Mask

SIZE	COLOUR	CODE
S	Beige	NOA-901-S-H
M	Beige	NOA-901-M-H
L	Beige	NOA-901-L-H
XL	Beige	NOA-901-XL-H

NOA-003

Seamless Cup Basic Bra

SIZE	COLOUR	CODE
32 B/C	Black	NOA-003-32 B/C-B
	Beige	NOA-003-32 B/C-H
34 B/C	Black	NOA-003-34 B/C-B
	Beige	NOA-003-34 B/C-H
34 D/E	Black	NOA-003-34 D/E-B
	Beige	NOA-003-34 D/E-H
36 B/C	Black	NOA-003-36 B/C-B
	Beige	NOA-003-36 B/C-H
36 D/E	Black	NOA-003-36 D/E-B
	Beige	NOA-003-36 D/E-H
38 B/C	Beige	NOA-003-38 B/C-H
38 D/E	Black	NOA-003-38 D/E-B
	Beige	NOA-003-38 D/E-H
40 B/C	Black	NOA-003-40 B/C-B
	Beige	NOA-003-40 B/C-H
40 D/E	Black	NOA-003-40 D/E-B
	Beige	NOA-003-40 D/E-H
42 B/C	Black	NOA-003-42 B/C-B
	Beige	NOA-003-42 B/C-H
42 D/E	Black	NOA-003-42 D/E-B
	Beige	NOA-003-42 D/E-H

NOA-007

Seamless Cup Bra with Implant Stabiliser

SIZE	COLOUR	CODE
32 B/C	Beige	NOA-007-32 B/C-H
34 B/C	Beige	NOA-007-34 B/C-H
34 D/E	Black	NOA-007-34 D/E-B
	Beige	NOA-007-34 D/E-H
36 B/C	Beige	NOA-007-36 B/C-H
36 D/E	Black	NOA-007-36 D/E-B
	Beige	NOA-007-36 D/E-H
38 B/C	Black	NOA-007-38 B/C-B
	Beige	NOA-007-38 B/C-H
38 D/E	Black	NOA-007-38 D/E-B
40 B/C	Black	NOA-007-40 B/C-B
	Beige	NOA-007-40 B/C-H
40 D/E	Black	NOA-007-40 D/E-B
42 B/C	Black	NOA-007-42 B/C-B
	Beige	NOA-007-42 B/C-H
42 D/E	Black	NOA-007-42 D/E-B

NOA-510

Implant Stabiliser Band

SIZE	COLOUR	CODE
ONE SIZE	Black	NOA-510-U-B
ONE SIZE	White	NOA-510-U-W

NOA-310

Women's Vest with Open Bust

SIZE	COLOUR	CODE
2XS	Black	NOA-310-2XS-B
	Beige	NOA-310-2XS-H
XS	Black	NOA-310-XS-B
	Beige	NOA-310-XS-H
S	Black	NOA-310-S-B
	Beige	NOA-310-S-H
M	Black	NOA-310-M-B
	Beige	NOA-310-M-H
L	Black	NOA-310-L-B
	Beige	NOA-310-L-H
XL	Black	NOA-310-XL-B
	Beige	NOA-310-XL-H
2XL	Black	NOA-310-2XL-B
	Beige	NOA-310-2XL-H
3XL	Black	NOA-310-3XL-B
	Beige	NOA-310-3XL-H
4XL	Black	NOA-310-4XL-B
	Beige	NOA-310-4XL-H

NOA-147

Women's Girdle with High-Back Calf Length

SIZE	COLOUR	CODE
S	Black	NOA-147-S-B
	Beige	NOA-147-S-H
M	Black	NOA-147-M-B
	Beige	NOA-147-M-H
L	Black	NOA-147-L-B
	Beige	NOA-147-L-H
XL	Black	NOA-147-XL-B
	Beige	NOA-147-XL-H
1XL	Black	NOA-147-1XL-B
	Beige	NOA-147-1XL-H
2XL	Black	NOA-147-2XL-B
	Beige	NOA-147-2XL-H
3XL	Black	NOA-147-3XL-B
	Beige	NOA-147-3XL-H

NOA-157

Women's Girdle with High-Back Short Length

SIZE	COLOUR	CODE
2XS	Black	NOA-157-2XS-B
	Beige	NOA-157-2XS-H
XS	Black	NOA-157-XS-B
	Beige	NOA-157-XS-H
S	Black	NOA-157-S-B
	Beige	NOA-157-S-H
M	Black	NOA-157-M-B
	Beige	NOA-157-M-H
L	Black	NOA-157-L-B
	Beige	NOA-157-L-H
XL	Black	NOA-157-XL-B
	Beige	NOA-157-XL-H
1XL	Black	NOA-157-1XL-B
	Beige	NOA-157-1XL-H
2XL	Black	NOA-157-2XL-B
	Beige	NOA-157-2XL-H
3XL	Black	NOA-157-3XL-B
	Beige	NOA-157-3XL-H
4XL	Black	NOA-157-4XL-B



For more information about our products please visit:
www.gcaesthetics.com or email info@gcaesthetics.com

Global Consolidated Aesthetics Limited, registered in Ireland with number 450181.
Registered office: Suite 601, Q House, Furze Road, Sandyford Industrial Estate,
Dublin 18, Ireland. © Copyright 2019 GC Aesthetics. A50-193.

REFERENCES

1. Gladfelter, J (2007). Compression Garments 101. Plastic Surgical Nursing: Vol. 27, Issue 2.
2. Illouz, Y.G (1989). Body sculpturing by lipoplasty. New York, NY: Churchill Livingstone.
3. Marena Data on File.
4. Fortes et. al (2008). Post Liposuction Compression Garment Comparison: Is there a difference?
A clinical study by the Aesthetics Center for Plastic Surgery.